



## Monoclonal Antibody to Fibrinogen

RMA03100 100ul

### **[ PROPERTIES ]**

**Source:** Monoclonal antibody preparation

**Host:** Rabbit

**Purification:** Protein A and Protein G affinity chromatography.

**Traits:** Liquid

Concentration: 0.5mg/ml

Ig Isotype: IgG

**UOM:** 100µl

**Cross Reactivity:** Human, Mouse, Rat

**Applications:** WB, IHC

### **[ IMMUNOGEN ]**

Immunogen: Recombinant protein of Fibrinogen.

### **[ APPLICATIONS ]**

WB: 1:1000

IHC: 1:1000

Optimal working dilutions must be determined by end user.

### **[ FORMULATION ]**

**Form & Buffer:** Supplied as solution form in 0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

### **[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 4°C for frequent use.

Aliquot and store at -20°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was

determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

**[ IDENTIFICATION ]**




Figure. Western Blot; Sample:  
Recombinant Fibrinogen, Human.

**[ IMPORTANT NOTE ]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.